


PLATE SEAL

ChemWell® Factory Installed Option



*Compatible with ChemWell®
2900 series analyzers*

Designed as an optional factory installed enhancement for use with specific ChemWell® brand analyzers, manufactured by:



It is recommended for automatic ELISA processing when the manual test instructions indicate sealing or covering the plate during one or more of the incubation steps. Since the automatic processor cannot peel plastic seals on and off, this novel device adapts to the analyzer and automatically seals and unseals the plate with no user intervention.

Many assays are not sensitive to whether or not they are covered but when the seal is required to match manual results,



Eliminates the need for films or tape